

DNA TESTING WITH ANTAGENE



EN

SAMPLING KIT

FOR OFFICIAL RECOGNITION THE SAMPLE MUST BE AUTHENTICATED BY A VETERINARIAN

To get an official recognition of the results (kennel club, breed registers) To certify the identity of the sampled animal

HOW TO TAKE A CHEEK SWAB

SAMPLING The cheek swab is a sampling of mucosal fragments

- Slide brush between the animal's cheek and gum
- Collect the sample by rolling and rubbing firmly the cheek swab on the inside of the pet's cheek for 20 seconds



Cheek

swab

PACKAGING

- Place the swab within the ethanol tube, if ethanol is missing, complete with modified alcohol
- Cut the swab just above the brush
 - Close the tube and make sure it does not leak
- Contact us to get a free sampling kit
- Copy the tube's number on the sampling certificate
- Check that all animal's information are accurate, then sign and stamp the sampling certificate
- Place the sample in the enveloppe provided for this purpose

Caution

To avoid contamination that would lead to the failure of the test, the animal must not eat, drink, suckle, nor lick anything or any other animal for 30 minutes prior taking the sample.

When no sampling kit is available

- Collect 1mL to 2mL of whole blood on EDTA
- → Label the tube with the animal name or its identification number (microchip or tattoo)
- Place the sample within an hermetic bag, suitable for sending a blood sample

HOW TO ORDER A SAMPLING KIT

Contact ANTAGENE



By email contact@antagene.com



Online

antagene.com/en/my-space



By phone +33 4 37 49 90 03 PAYMENT





Online

antagene.com/en/my-space

- · When you order your sampling kit
- Or on invoice receipt, once your sample is received by ANTAGENE

RESULTS

PDF document sent by email



To the owner and veterinarian

PDF document uploaded online



Online

antagene.com/en/my-space



SAMPLING CERTIFICATE

Tube No.

Date :	For ANTAGENE use or	
Position:		Version :

	☐ Blood sample (EDTA)	☐ Other :
Conv the tube number shove	Identify the tube (name or chin)	

	SAMPLE AUTHENTICATION BY A VETERINAR Authentication can not be performed re			OFFIC	CIAL	RES	ULT				
	I, the undersigned Dr. LAST NAME First Clinic Name Address										
	ZIP Code City & State Clinic email for results Invoicing email					Country					
ERINARIAN	□ Dog □ Cat Breed Registered Name Usual Name □ I certify having sampled this animal myself, after having checker										
VETER	☐ By microchip scanning ☐ By tattoo reading Identification Number										
31/12/2025	REQUESTED ANALYSIS □ 45€ DNA profile (ISAG 2006 / ISAG F2014) Doesn't include the Cat Genetic Blood Group										
Faxes included (VAT), valid until 31/12/2025	 ☐ 75€ One Test Specify the test ☐ 105€ DNA profile + One Test Specify the test ☐ 168€ Dog Genetic Check-Up or Cat Genetic Check-Up 	□ 75€ MDR1 □ 108€ CEA □ 108€ prcd-F □ 120€ SH if older than 3 months blood sample mandatory Included in the Dog Genetic Check-Up □ 75€ Cat Genetic Blood Group (all breeds)								- 11	
Тахе	Dog Genetic Check-Up or Cat Genetic Check-Up + DNA profile Veterinarian's sta Sampling date D D M M Y Y The sample must be sent under the the veterinarian's liability.	Inc	luded in the	The	Vete	rinaria ian accep	ots the c	urrent 1	Terms and		
ADMINISTRATION	INVOICE TO Owner (email required) Veterinarian (email required) Other: LAST NAME, First name, phone number, email, postal address No results will be delivered in the absence of payment. Results are automatically sent by emails	to the anir	nal's owner a		☐ Alro	k card eady pa ocess yment u	iid durir Ipon inv	ng the voice r	·		
OWNER	☐ Mr. ☐ Mrs. LAST NAME E-mail Address ZIP Code City & State Country Phone number ☐ I confirm that informations concerning the animal and the requested analysis			The ow	ner certi	wner's	s signa he/she o	iture owns th Condit	ne animal a ions of sale		
	☐ I agree that the result can not be recognized as official in the absence of anim			authent	ication	by a v	eterin	arian			